

### Chlorophyll Fluorometer Characterization in Uranine liquid Proxy

**Date:** 12/22/15  
**Serial #:** AFLT-015  
**Tech:** SML

**Dark Counts** 0.2091 volts  
**CEV** 1.568volts  
**SF** 15.58

**FSV** 5.36 volts

**Linearity:** 0.999 R<sup>2</sup> (0–1.5 volts)  
0.995 R<sup>2</sup> (0– 5.45 volts)

Notes:

**Dark Counts:** Signal output of the meter in clean water with black tape over detector.

**CEV** is the chlorophyll equivalent voltage. This value is the signal output of the fluorometer when using a Uranine dye fluorescent proxy that has been determined to be approximately equivalent to **21.17 µg/l** of a *Thalassiosira weissflogii* phytoplankton culture.

**SF** is the scale factor used to derive chlorophyll concentration from the signal voltage output of the fluorometer. The scale factor is determined by using the following equation:  $SF = (21.17) / (CEV - \text{dark})$ .

**FSV** is the maximum signal voltage output that the fluorometer is capable of.

Chlorophyll concentration expressed in µg/l ( $\text{mg/m}^3$ ) can be derived by using the following equation: ( $\mu\text{g/l}$ ) =  $(V_{\text{measured}} - \text{dark}) * SF$

The relationship between fluorescence and chlorophyll-*a* concentrations in-situ is high variable. The scale factor listed on this document was determined by using a mono-culture of phytoplankton (*Thalassiosira weissflogii*). The population was assumed to be reasonably healthy and the concentration was determined by using the absorption method. To accurately determine chlorophyll concentration using a fluorometer you must perform secondary measurements on the populations of interest. This is typically done using extraction based measurement techniques on discrete samples. For additional information on determination of chlorophyll concentration see [ Standard Methods For The Examination Of Water And Wastewater] part 10200 H published jointly by: American Public Health Association, American Water Works Association and Water Environment Federation.

PO Box 518  
620 Applegate St.  
Philomath OR 97370  
wetlabs@wetlabs.com



(541) 929-5650  
Fax (541) 929-5277  
www.wetlabs.com

**Date: 12/22/15**

**Customer: University of Alaska/ Steve Hartz**

**S/N# AFLT-015**

**Technician: SML**

**Repairs and Modifications: Evaluated instrument. Found instrument to be responsive and current draw to be normal. Tested then characterized instrument in Uranine proxy. Replaced case seals.**

**Comments: New Characterization sheet included.**